# Anclacomp nano

## Instrucciones Composite Anclacomp Nano

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# Anclacomp nano

## Instructions for use

Nano Composite is a light-curing, high-lustre-polish hybrid composite with an ultrafine, radiopaque porcelain filler for use in adhesive filling treatment. It can be polished to a high lustre.

Due to the ultra-fine particle filler, extremely homogeneous restorations can be placed which are easily polished to a high lustre. The chameleon effect matches the shade of the filling perfectly to the tooth structure.

Nano Composite is supplied in practical screw syringes.

## Composition:

**Monomer matrix:** diurethane dimethacrylate, butanediol dimethacrylate, isopropylidene-bis[2(3)-hydroxy-3(2)-(4-phenoxy)propyl] bismethacrylate [Bis-GMA]

**Total filler content 83% by weight (70% by volume)**

- Glass filler (medium grit size 0.7 µm)
- Pyrogenic silica (medium grit size 12 µm)
- Agglomerated nanoparticles (medium grit size 0.6 µm)

## Indications

- Direct anterior and posterior restorations in Black's class I, II, III, IV and V cavities.
- Indirect restorations such as inlays, onlays and laminate veneers
- Extended fissure sealing in molars and premolars
- Cores
- Splinting mobile teeth
- Adjusting the contours and shade to improve aesthetics

## Application

### Preparatory measures

Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita® shade guide to select the shade while the tooth is still moist.

- Cavity preparation**  
Prepare the cavity minimally invasively as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.
- Pulp protection / Cavity liner**  
If an enamel-dentine adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.
- Approximal contact areas**  
When filling cavities with approximal sections, place a transparent matrix and fix it in place.
- Etching with Etching**  
Apply Etching first to the enamel areas of the cavity and let it take effect for 30 seconds. Then fill the whole cavity (dentine) with etching gel and let it take effect for a further 15 seconds. The etching time in the dentine should not exceed 20 seconds. Then rinse the cavity thoroughly with a water spray and dry it with oilfree compressed air. Do not overdry the dentine. Dried, etched enamel surfaces have a chalky-white appearance and must not be contaminated before the bonder is applied. If the surface becomes contaminated with saliva, rinse and dry again and re-etch if necessary.
- Application of the bonding agent**  
Refer to the Anclabond instructions for details.
- Placing the composite in the cavity**  
Remove the required amount of composite from the screw syringe, insert it into the cavity using standard metal instruments and contour. Layers should not be more than 2 mm thick. A dispersion layer,

a thin, non-polymerised film, remains on the surface after curing due to the effect of the oxygen in the air. This produces the chemical bond between the layers and must not be touched or contaminated with moisture. Curing

The exposure time for all shades per layer is 40 seconds with a standard halogen polymerisation unit or an LED polymerisation lamp or 2 x 3 seconds with a plasma polymerisation unit. Hold the light guide as close as possible to the surface of the filling. Polymerise each surface of complex fillings.

### 7. Trimming

Nano Composite can be trimmed and polished immediately after curing using finishing diamonds, flexible disks, silicone polishers and polishing brushes. Check the occlusion and articulation and spot grind to eliminate high spots or undesirable paths of articulation from the surface of the filling.

## Indirect method Cavity preparation

The cavity should be prepared as minimally invasively as possible with only slightly diverging sides. All internal line and point angles must be rounded. Avoid slice preparations. Prepare a flat cervical shoulder – do not bevel it. Any unavoidable undercuts must be blocked out with glass ionomer cement. Use slightly tapering diamonds with rounded tips for the preparation. Coat those areas of dentine in close proximity to the pulp with a thin layer of calcium hydroxide material. Cavity liners containing eugenol are contraindicated.

## Impression and temporary restoration

Once the impression has been taken e.g. with Anclasil Oneputty, a composite temporary restoration is fabricated e.g. with Anclaprov Tempo. This must only be cemented with a non-eugenol cement.

## Fabricating an inlay

Cast the impression with high strength dental stone in the laboratory. Allow the model to set and pull off the impression. Block out the under-cuts and apply an oil-free separating agent to the model. Build up the inlay layer-by-layer on the model. Build up the approximal and deep occlusal sections first. Each individual layer must not be more than 2 mm thick and is cured separately with a commercially available light curing lamp. The finished inlay is then released from the die and cured fully with polymerisation standad lamp. Trim the occlusal surface with fissure burs and polish to a high lustre with silicone polishers and diamond paste. Clean the inlay thoroughly with soap and water, spray with air / water to rinse and dry.

## Fitting inlays, onlays or veneers

Remove the temporary restoration and clean the cavity. Place a rubber dam, then clean and dry the prepared surface of the tooth. Check the fit of the restoration applying light pressure. Do not force the restoration to fit. Grind the fitting surface if necessary to improve the fit. The occlusion and articulation should not be checked during the try-in of the inlay, as there is a risk of fracturing the inlay. Etch the tooth with Anclatetch (refer to the Anclatetch instructions for use), rinse thoroughly with the water spray and dry the etched surfaces with oil-free compressed air. The dried, etched enamel surface has a chalky-white appearance and must not be contaminated before applying the bonding Anclabond. If the surface is contaminated with saliva, rinse and dry again and re-etch if necessary. Apply a thin layer of Anclabond bonding agent with a brush to the etched enamel and cavity walls, rub it in and cure for 40 seconds with a standard halogen polymerisation unit. After curing, a dispersion layer forms that should not be removed, as it produces the chemical bond to the filling material. Sandblast the interior of the composite-restoration surface and the clean with alcohol, then apply Anclabond without light curing. Mix commercial dual-hardening composite cement and apply to the interior surface of the restoration. Press the inlay carefully into position and remove larger amounts of excess with a spatula. Remove approximal excess with a probe or dental floss. Press lightly on the inlay with a ball-end stopper until curing is complete to ensure that the inlay remains in the correct position. Beginning at the approximal, cure each surface of the composite for 40 seconds with a standard halogen polymerisation unit or LED polymerisation unit. Remove excess with fine-grit diamond rotary instruments and then diamond-coated finishing strips.

Check the occlusion and adjust it if necessary. Complete final finishing and polishing with a finishing and polishing set.

**Note: If an inlay is more than 2 mm thick, a dual-curing composite must be used.**

## Please note

- The working time beneath a dental light is approximately 2 minutes.
- When placing time consuming restorations, to prevent the composite curing prematurely the dental light should be moved away from the site temporarily or the composite covered with foil impervious to light.
- Use a light curing unit with an emission spectrum of 350 – 500 nm for curing this material. As the required physical properties can only be achieved if the lamp is functioning correctly, its luminous intensity must be checked regularly as described by the manufacturer.

## Storage

Do not store above 25°C. Avoid direct sunlight. Close the screw syringes tightly immediately after use. The material should be at room temperature before use. Retract the plunger of the syringe slightly to prevent the apertures becoming blocked. Do not use after the expiry date (refer to label on syringe). For use by dentists only. Keep out of reach of children. This product was developed specifically for the described range of applications. It must be used as described in the instructions. The manufacturer is not liable for damage caused by handling or processing the material incorrectly.

## Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. To prevent possible reactions of the pulp in cavities where the dentine is exposed, the pulp must be protected adequately (e.g. calcium hydroxide preparation).

## Contra-indications / interactions

If a patient has known allergiens against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product. Unpolymerized composite may cause skin allergies. The user must take adequate precautions.

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners (such as zinc-oxide eugenol cements) containing such substances.

\* Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

## Trouble shooting

Problem	Cause	Remedy
Composite does not cure completely	Light output of the light-curing lamp inadequate	Check the light output and change the light source if required
	Emitted wavelength range of the lightcuring lamp inadequate	Consult the manufacturer of the lightcuring lamp. Recommended wavelength range: 350 - 500 nm
Composite in the syringe is sticky and soft, colourless liquid separates in the syringe	Material has been stored over a longer period at ≥ 25°C	Adhere to storage temperature; store in a refrigerator
	Material has been kept too long in a syringe warmer	Never keep a syringe longer than one hour per application in a syringe warmer
Composite appears too hard and firm in the syringe	The material has not been heated to room temperature after removal from the refrigerator	Allow the composite to heat to room temperature before use; use a syringe warmer if necessary
	Syringe not properly sealed, composite partially cured	Always seal the syringe properly with the cap after removing composite
Inlay / onlay is not properly retained when fitted	Restoration is too opaque to be cemented using only lightcuring composite	Use dual-curing luting composite
Composite does not cure all the way through (dark or opaque shades)	Composite layers applied too thickly for each curing cycle	Adhere to a max. layer thickness of 2.0 mm per increment
Restoration appears too yellow compared with the shade guide	Inadequate curing of the composite layer	Repeat the exposure cycle several times; min. 40 sec.

Ancladén, S.L.  
Rubinstein, 2  
08022 Barcelona (España)  
www.ancladen.com

0434