# Anclacomp bulkfill

## Instrucciones Composite Anclacomp-Bulkfill

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# Anclacomp bulkfill

## Instructions *Bulk-Fill Composite*

Bulk-Fill Composite is a light curing, flowable, highly radiopaque (210 % Al), composite for posterior restorations, for direct filling and the "bulk fill" technique.

### Composition

- Monomer matrix: diurethane dimethacrylate, butanediol dimethacrylate Total
- Filler: 77 % (by weight) glass filler, pyrogenic silicic acid

### Indications

- Fillings with layer thicknesses up to 4 mm in Black's class I and II cavities
- Cavity lining - as a first (bottom) layer in Black's class I and III cavities

### Contraindications

Unpolymerized composite may cause skin allergies. The user must take adequate precautions (e.g. protective gloves). Do not use this material, if the user or patient should have known allergies against one of the constituents listed under "Composition".

### Side-effects

To prevent possible reactions of the pulp in cavities where the dentine is exposed, the pulp must be protected adequately (apply e.g. a calcium hydroxide preparation).

### Interaction with other substances

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners containing such substances.

### Application Preparatory measures

Before commencing the treatment, clean the tooth with non-fluoride polishing paste.

Prepare the cavity minimally invasively as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.

### Approximal contact areas

For cavities with approximal sections, place a transparent matrix and fix it in place.

### Etching with Etching Gel

Ensure a very precise application of Etching Gel, beginning with the enamel margins. Allow it to react for at least 15 seconds before spraying thoroughly with water and drying the etched surfaces with oil-free compressed air. Avoid dehydrating the dentine. After drying, the etched surfaces of the enamel appear chalky/white and must not be contaminated prior to applying the bonding agent. Should these surfaces be

contaminated with saliva, rinse and dry them again – re-etch if necessary.

### Applying the bonding agent (Bond LC)

Refer to the Bond LC instructions for details.

### Application of Bulk-Fill Composite

Place layers of max. 4 mm of Bulk-Fill Composite directly into the cavity using the curved application tips supplied with the material. For hygiene reasons, these application tips are must only be used once and then discarded. When applying the material with the tip, make certain that no air bubbles become entrapped in the material. Wet the prepared tooth surfaces thoroughly. Polymerize each layer with a conventional polymerization unit (e.g. Translux CL) for 40 seconds, holding the light guide as close as possible to the surface. During polymerization, a dispersion layer forms on the material surface. This layer must not be touched or removed. It forms the bond with the layer of composite applied subsequently (e.g. with a universal or posterior composite) or for creating aesthetic cusps. Alternatively, Bulk-Fill Composite can be applied in a thin layer to form a liner on exposed dentin. The material can be polymerized in layers of up to 4 mm thickness.

### Finishing

Bulk-Fill Composite is finished and polished with finishing diamond tools, flexible discs, silicone polishers and polishing brushes immediately after polymerization.

### Please note

When creating large and time-consuming restorations either cover the material with a lightproof foil or temporarily move the operating light away from the site in order to avoid premature polymerization of the material.

For hygienic reasons, the curved application tips supplied with the material must be used only once!

Use a polymerization unit with an emission spectrum of 350–500 nm to cure this material. As the required physical properties can only be achieved if the lamp is functioning correctly, its luminous intensity must be checked regularly as described by the manufacturer.

### Storage

Store at temperatures between 10-25°C. Avoid direct sunlight. Close the syringes tightly immediately after use. Let the material reach room temperature before use. After use, retract the plunger of the syringe slightly to prevent the apertures from becoming blocked. Do not use after the expiry date (refer to label on syringe). For use by dentists only. Keep out of children's reach. This product was developed specifically for the described range of application. It must only be used in accordance with the above instructions. The manufacturer will not be held liable for damages caused by incorrect handling or incorrect processing of the material.

### Trouble shooting

Problem	Cause	Remedy
Composite does not cure completely	Light output of the light-curing lamp inadequate	Check the light output and change the light source if required
	Emitted wavelength range of the lightcuring lamp inadequate	Consult the manufacturer of the lightcuring lamp. Recommended wavelength range: 350 - 500 nm
Composite in the syringe is sticky and soft, colourless liquid separates in the syringe	Material has been stored over a longer period at $\geq 25^{\circ}\text{C}$	Adhere to storage temperature; store in a refrigerator
	Material has been kept too long in a syringe warmer	Never keep a syringe longer than one hour per application in a syringe warmer
Composite appears too hard and firm in the syringe	The material has not been heated to room temperature after removal from the refrigerator	Allow the composite to heat to room temperature before use; use a syringe warmer if necessary
	Syringe not properly sealed, composite partially cured	Always seal the syringe properly with the cap after removing composite
Inlay / onlay is not properly retained when fitted	Restoration is too opaque to be cemented using only lightcuring composite	Use dual-curing luting composite
Composite does not cure all the way through (dark or opaque shades)	Composite layers applied too thickly for each curing cycle	Adhere to a max. layer thickness of 2.0 mm per increment
Restoration appears too yellow compared with the shade guide	Inadequate curing of the composite layer	Repeat the exposure cycle several times; min. 40 sec.

 **Ancladén**  
Ancladén, S.L.  
Rubinstein, 2  
08022 Barcelona (España)  
www.ancladen.com