# Anclacem implant

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E

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Editing date: 13.09.2010

Instructions for use

## Resin-based temporary luting cement for implant-retained restorations

Dual-curing semi-permanent implant cement with very low displacement resistance and thus extremely low film thickness. Due to highly cross-linked, ductile polymer structure, high compression strength and hermetic sealing are achieved. Excess material can easily be removed in large pieces (no crumbling) during plastic phase. Luted restorations are effortlessly removed without damage with standard devices (Crown Butler, Corona Flex). Compatible with all conventional material combinations, i.e. zirconium and metal frameworks, ceramics, titanium. Antibacterial, eugenol-free, odourless and tasteless; natural opaque colour for screening metallic abutments.


Fig. 1

### 1. Preparation and cleaning of implant abutments and inner surfaces of restorations

Carefully check fit, occlusion and proximal contact points. Complete all finishing procedures. Degrease, clean (e.g. alcoholic solutions), dry and isolate prepared abutments. Make sure the interior surfaces of the superstructure are clean and dry (Fig. 1).


Fig. 2

### 2. Mixing and applying

Bleed the syringe before installation of the automix cannula until material emerges uniformly from both openings (Fig. 2). Insert the plunger in the body of the cartridge. Remove the cap from the cartridge by turning clockwise and attach mixing cannula. Pay attention to the guides on cannula and cartridge. Lock in place by turning clockwise (Fig. 3). Again bleed small amount to ensure uniform and even material flow. Due to finely tuned dispensing pressure, dosing can be performed individually and precisely. Dispense the material with uniform pressure. After use, leave mixing cannula on cartridge until next application to keep cement from drying out.

Fig. 3

### 3. Inserting the restoration

Apply a thin layer of **Anclacem Implant** directly into the superstructure (Fig. 4). Carefully position the restoration on the implant abutment with slight pressure (Fig. 5); then press down firmly. Remove excess material during the gel phase, after around 2 – 3 minutes, with a suitable instrument (Fig. 6).

Fig. 4

### 4. Working time / light curing

The working time at room temperature is around 80 seconds. The cement is cured at body temperature in about 5 – 6 minutes. Direct exposure to light accelerates setting time. Light curing (marginal region, excess material) is achieved with approx. 20 seconds per surface (Fig. 7). When **Anclacem Implant** is used with translucent materials (e.g. thin layers of ceramic provisional crowns and bridges), increase exposure to light through the material to 30 - 40 seconds per surface.


Fig. 5

### Precautions

- Do not damage the underlying titanium surface when removing excess material. If necessary, use plastic instruments (implant scalers).
- Leave no cement residues in the mouth or gingival sulcus.
- Use no adhesion-reducing additives. They could uncontrollably render the retentive qualities of **Anclacem Implant**.
- Do not manually remove or mix cement components.
- Avoid contact of unmixed components with skin or mucosa.
- In very rare cases, **Anclacem Implant** may cause sensitization. Should such reaction occur, discontinue use of this material.
- Observe the safety data sheet.
- Do not use **Anclacem Implant** for the luting of hybrid restorations with natural teeth and implant abutments.
- Only for intended use by qualified personnel.
- Do not rough or sand-blast abutment surface.


Fig. 6

### Autopolymerisation:



Fig. 7

### Light-curing:

### Indication for Use:

- Temporary luting of implant-retained restorations

### Technical data:

- Mixed volume:** 5 ml (mini-mix)
- Mixing ratio:** 4:1
- Colour code:** Base: white-opaque  
Catalyst: semi-transparent
- Mixing time:** None: automix (mini-mix)
- Working time:** approx. 80 sec.
- Gel phase:** 2 - 3 min. \* (for removal of excess material during placement)
- Setting time in the mouth:** 5 - 6 min. \*
- Setting time with light-curing:** approx. 20 sec. per surface; 30-40 sec. through translucent materials
- Film thickness:** < 10 µm
- Compression strength:** approx. 85 MPa
- Application temperature:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- Storage:** below 25° C / 77° F

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Higher temperatures accelerate, lower temperatures retard a. m. times.

### Ordering information:

**Anclacem Implant Standard packing A-050001**  
cartridge of 5 ml mini-mix 4:1  
10 mixing cannulas

Mixing cannulas 4:1  
25 pcs. **A-050002**

# Ancladén

Ancladén, S.L. | Rubinstein, 2 Entlo 1ª | 08022 Barcelona  
93 417 29 84 | ancladen@ancladen.com | www.ancladen.com